

**OPIA**  
INNOVATION



 **OPIA**  
*Pirmax*

## Standard Electric Wheelchairs



PA200



PA201



PA202



PA204



PA204K



PA207



PA208



PA210



PA303

## Manual Wheelchairs



PM110



PM111



PM112



PM210



PM220



PM220A



PM230



PM240



PM250



PM320



PM321



PM330



PM331



PM500



PM930



PM930L



PM931



PM931L



PM9070



PM9031C




BUGGY



“

**Standard  
Electric Wheelchair**



Chassis	Foldable with Battery Removed
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	100-120 kg
Wheel Dimensions	8 "x2" (Front) 12.5 "x2.25" (Rear)
Brake System	Electromagnetic Brake
Batteries	Gel 12V33AH*2 / 32A-33A
Charger	DC24V 4A - (5A/6A)
Charging Time	6-8h
Motor	24V/250W*2
Joystick	24V/40A / 50A
Weight without Battery	49 kg
Model	Suitable for 18" 20" and 22" chassismodel construction.



62-63 cm




44-46 cm



100-120 kg



60-63 kg



32/33



250 W



40/50



Chassis	Foldable with Battery Removed
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	120-135 kg
Wheel Dimensions	8 "x2" (Front) 12.5 "x2.25" (Rear)
Brake System	Electromagnetic Brake
Batteries	Gel 12V32AH*2-33A
Charger	DC24V5A (5A/6A)
Charging Time	6-8h
Motor	24V/320W*2 - 340W*2
Joystick	24V/50A - 50A/70A
Weight without Battery	59 kg
Model	Suitable for 18" 20" and 22" chassismodel construction.



62-63 cm




44-50 cm



120-135 kg



57-59 kg



32/35 Ah










340 W



50/70



Chassis	Foldable with Battery Removed
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	120-135 kg
Wheel Dimensions	8"/10" (Front) - 12.5"/14" (Rear)
Brake System	Electromagnetic Brake
Batteries	Gel 12V33AH*2 / 40AH
Charger	DC24V5A (5A/6A)
Charging Time	6-8h
Motor	24V/320W*2 - 340W*2
Joystick	24V/50A - 70A
Weight without Battery	60 kg
Model	Suitable for 18" 20" and 22" chassismodel construction.

						
60-63 cm	44-54 cm	120-135 kg	60-62 kg	5/6 Ah	33-40 Ah	320/340 W



Chassis	Aluminum Chassis - Folding
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	120-135
Wheel Dimensions	8 "x2" (Front) 12.5 "x2.25" (Rear)
Brake System	Electromagnetic Brake
Batteries	Gel 12V32AH*2 - 33AH
Charger	DC24V4A (5A/6A).
Charging Time	6-8h
Motor	24V/320W*2 or 24V/340W*2
Joystick	24V/50A - 70A o
Weight without Battery	65 kg
Model	Suitable for 18" 20" and 22" chassismodel construction.



45-50 cm

45-50 cm

120/135kg

52-55 kg

32/33 Ah

320/340W



Chassis	Aluminum Chassis - Folding
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	100-120 kg
Wheel Dimensions	8 "x2" (Front) 12.5 "x2.25" (Rear)
Brake System	Electromagnetic Brake
Batteries	Gel 12V32AH*2 - 33AH
Charger	DC24V4A (5A/6A)
Charging Time	6-8h
Motor	24V/320W*2 - 24V/340W*2
Joystick	24V/50A - 70A
Weight without Battery	65 kg
Model	Suitable for 14" 16" and 18" chassis model construction.



50-62 cm




40-50 cm



100/120kg



44-49 kg



32/33 Ah



320/340W





Chassis	Aluminum Chassis - Folding
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	100-120 kg
Wheel Dimensions	8 "x2" (Front) 12.5 "x2.25" (Rear)
Brake System	Electromagnetic Brake
Batteries	Lithium 24V 20AH
Charger	DC24V3A (5A/6A)
Charging Time	5 h
Motor	24V/250W*2
Joystick	24V/50A
Weight without Battery	24 kg
Model	Suitable for 14" 16" and 18" chassismodel construction.



55-60 cm




44-50 cm



100-120 kg



20-25 kg



24 Ah



250 W



Chassis	Steel Chassis - Folding
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	120-135 kg
Wheel Dimensions	8"/10" (Front) - 12.5"/14" (Rear)
Brake System	Electromagnetic Brake
Batteries	Ge 12V32AH*2 - 33A/40A
Charger	DC24V4A (5A/6A)
Charging Time	6-8h
Motor	24V/320W*2 - 340*2
Joystick	24V/50A - 70A
Weight without Battery	65 kg
Model	Suitable for 18" 20" and 22" chassismodel construction.



55-65 cm




45-50 cm



120-135 kg



65/75 kg



32-33-40



320-340 W



Chassis	Steel - Shock Absorber
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	125-145 kg
Wheel Dimensions	8 "x2" (Front) 12.5 "x2.25" (Rear)
Brake System	Electromagnetic Brake
Batteries	Ge12V32AH*2 - 33
Charger	DC24V5A (5A/6A)
Charging Time	6-8h
Motor	24V/320W*2 - 340W*2
Joystick	24V/50A - 70A
Weight without Battery	72 kg
Model	Suitable for 18" and 20" chassismodel construction.



62-63 cm




44-50 cm



125-145 kg



60-65 kg



32-33 Ah



320/340 W



Chassis	Steel, Captain Seats
Maximum Speed	6-8 km/h
Max. Driving Angle	12°
Max. Driving Range	25±5km
Transport Capacity	150-165 kg
Wheel Dimensions	265x85 front - 300-8" rear
Brake System	Electromagnetic Brake
Batteries	Ge 12V40AH*2 - 56-64
Charger	DC24V5A (5A/6A)
Charging Time	10-12 h
Motor	24V/340W*2 - 450W
Joystick	24V/50A -70A
Weight without Battery	71kg
Model	Suitable for 18" 20" and 22" chassismodel construction.



44-54 cm




65-70 cm



150/165 kg



60-65 kg



40-56-64



340-450 W

MADE IN  
TURKEY



# Transport WHEELCHAIR



 **OPiV**  
*Pirmax*







Chassis	Chromed Steel Chassis, Folding
Front Wheel	6 "x1.25" padding
Rear Wheel	24" filler (Aluminum composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Fixed, PVC
Footrest	May come out
Fabric	Black Leather
Backrest	Fixed
Others	Calf Band, Reflector, Hand Brake, Seat Belt
Model	Suitable for 18" 20" and 22" chassismodel construction.



55-60 cm




40-46 cm



100-120 kg



14-17 kg





Different Rim Options





Chassis	Chromed Steel Chassis, Folding
Front Wheel	6 "x1.25" padding
Rear Wheel	24" infill (Aluminum composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Fixed, PVC
Footrest	May come out
Fabric	Black Leather
Backrest	Fixed
Others	Baldır Bandı, Reflektör, El Freni, Emniyet Kemerı
Model	Suitable for 18" 20" and 22" chassismodel construction.

				
55-60 cm	40-46 cm	100-120 kg	14-17 kg	Different Rim Options



Chassis	Steel Chassis, Folding
Front Wheel	8 "x1.25" padding
Rear Wheel	24" infill (Aluminum composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Fixed, PVC
Footrest	Removable, Side Opening, Folding and High Adjustable
Fabric	Black Leather
Backrest	Fixed
Others	Calf Band, Reflector, Hand Brake, Seat Belt
Model	Suitable for 18" 20" and 22" chassismodel construction.



55-62 cm




44-46 cm



120-135 kg



16-19 kg



Different Rim Options



Chassis	Aluminum, Folding
Front Wheel	8 "x1.25" padding
Rear Wheel	24" filler
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, Side Opening, Folding and High Adjustable
Fabric	Black trim (flame retardant)
Backrest	Folding
Others	Calf Support, Reflector, Hand Brake, Seat Belt
Model	Suitable for 18" 20" and 22" chassis model construction.



50-55 cm




36-44 cm



120 -135 kg



15-17 kg



Different Rim Options



Chassis	Steel, Folding
Front Wheel	8 "x1.25" padding
Rear Wheel	24" filler (Aluminum composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Foldable and Detachable, PVC
Footrest	Removable, Side Opening, Folding and High Adjustable
Fabric	Black trim (flame retardant)
Backrest	Fixed
Others	Calf Band, Reflector, Hand Brake, Seat Belt
Model	Suitable for 18" 20" and 22" chassis model construction.

60-62 cm	40-46 cm	120-135 kg	15-17 kg	Different Rim Options



Chassis	Steel, Folding
Front Wheel	200x50 PU filling
Rear Wheel	24" infill (Aluminum composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, Side Opening, Folding and High Adjustable
Fabric	Black trim (flame retardant)
Backrest	Folding
Others	Waist Calf Support, Reflector, Hand Brake, Seat Belt
Model	18" 20" ve 22" şasselere uygun model yapımı.



50-65 cm




44-60 cm



120-150 kg



15-17 kg



Different Rim Options



Chassis	Folding Steel
Front Wheel	200x50 PU filling caster
Rear Wheel	24" infill (Aluminum, Composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, Side Opening, Folding and High Adjustable
Fabric	Black trim (flame retardant)
Backrest	Folding
Others	Calf Band, Reflector, Hand Brake, Safety Belt, Accompanying Brake, Seat Cushion
Model	Suitable for 18" 20" and 22" chassis model construction.

50-55 cm	40-46 cm	120-135 kg	15-17 kg	Different Rim Options



Chassis	Steel, Folding
Front Wheel	200x50 PU filling
Rear Wheel	24" infill (Aluminum, Composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Detachable, Side Opening, Folding and High Adjustable
Fabric	Black trim (flame retardant)
Backrest	Folding
Others	Calf Support, Reflector, Hand Brake, Safety Belt, Accompanying Brake, Seat Cushion
Model	Suitable for 18" 20" and 22" chassis model construction.



**50-65 cm**




**40-46 cm**



**120-150 kg**



**15-17 kg**



**Different Rim Options**



Chassis	Aluminum Folding
Front Wheel	8 "x2" padding
Rear Wheel	24" fill
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, side opening, folding and high adjustable
Fabric	Black trim (flame retardant)
Backrest	Fixed
Others	Calf band, reflector, hand brake, seat belt, companion brake, seat cushion
Model	Suitable for 18" 20" and 22" chassismodel construction.



60-65 cm



44-55 kg



120-135 kg



14-18kg





Chassis	Aluminum, Folding
Front Wheel	6 "x1.25" Padding
Rear Wheel	16" Filling
Seating Width	18" (46 cm)
Carrying Capacity	100-120 kg
Armrest	Height adjustable and removable.PU
Footrest	Detachable, folding and high adjustable and angle adjustable
Fabric	Black trim (flame retardant)
Backrest	Angle adjustable, Tilting
Others	Calf support, reflector, hand brake, seat belt, head, front and side support
Model	Suitable for 18" 20" and 22" chassismodel construction.



55-60 cm



44-50 cm



100-120 kg



17-20 kg




Chassis	Steel
Front Wheel	6 "x1.25" Padding
Rear Wheel	24" Filling
Seating Width	18" (46 cm)
Carrying Capacity	100-135 kg
Armrest	Height adjustable and removable, PU
Footrest	Removable, folding and high-set and angle adjustable
Fabric	Black trim (flame retardant)
Backrest	Angle Adjustable
Others	Calf support, reflector, hand brake, x type seat belt, head support, front and side support
Model	Suitable for 18" 20" and 22" chassismodel construction.

60-63 cm	46-50 cm	100-135 kg	17-22 kg	Different Rim Options



“  
**Pediatric  
Manual Wheelchair**”



Chassis	Steel, collapsible
Front Wheel	6 "x1.25" Padding
Rear Wheel	22" Filling
Seating Width	14" (36 cm)
Carrying Capacity	100-120 kg
Armrest	Fixed, PVC
Footrest	Fixed
Fabric	Plaid Flooring (flame retardant)
Backrest	Fixed
Others	Calf Band, Reflector, Hand Brake, Seat Belt
Model	Model construction suitable for 14" and 16" chassis.



50-55 cm



35-40 cm







100-120 kg



17-19 kg



Chassis	Aluminum Folding
Front Wheel	8 "x1.25" Padding
Rear Wheel	22" Filling
Seating Width	14" (36 cm)
Carrying Capacity	75-110 kg
Armrest	Can be folded back and out, PU
Footrest	Removable, side opening, folding and high adjustable
Fabric	Black trim (flame retardant)
Backrest	Folding
Others	Calf Band, Reflector, Hand Brake, Safety Belt, Accompanying Brake, Seat Cushion
Model	Model construction suitable for 14" and 16" chassis.

			
50-55 cm	30-40 cm	75-110 kg	13-14 kg







Chassis	Aluminum, Folding
Front Wheel	6 "x1.25" Padding
Rear Wheel	16" Filling
Seating Width	14" (36 cm)
Carrying Capacity	100-120 kg
Armrest	Height adjustable and removable,PU
Footrest	Detachable, folding and high adjustable and angle adjustable
Fabric	Black trim (flame retardant)
Backrest	Angle adjustable, Tilt
Others	Calf support, reflector, hand brake, seat belt, head, front and side support
Model	Suitable for 14" and 16" chassismodel construction.



55-60 cm



30-40 cm



90-110 kg




17-20 kg





Chassis	Aluminum, Folding
Front Wheel	6 "x1.25" padding
Rear Wheel	22" Filling
Seating Width	16" (40 cm)
Carrying Capacity	120-135 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, Side Opening, Folding and High Adjustable
Fabric	Black trim (flame retardant)
Backrest	Folding
Others	Calf Support, Reflector, Hand Brake, Seat Belt
Model	Suitable for 14" and 16" chassis model construction.



50-55cm

36-44 cm

120-135 kg

15-17 kg

Different Rim Options



Chassis	Steel, Folding
Front Wheel	6 "x1.25" padding
Rear Wheel	22" Filling
Seating Width	14" (36 cm)
Carrying Capacity	90-110 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, side opening, folding and high adjustable
Fabric	Black trim (flame retardant)
Backrest	Fixed
Others	Calf Band, Reflector, Hand Brake, Seat Belt
Model	Suitable for 14" and 16" chassis model construction.



50-55 cm



34-42 cm



90-110 kg



15-17 kg



Chassis	Aluminum, Folding
Front Wheel	6 "x2" Padding
Rear Wheel	10" Filler
Seating Width	16" (34cm)
Carrying Capacity	50 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, Angle Adjustable, Folding and High Adjustable
Fabric	Blue - Orange
Backrest	Can lie down, can rotate bilaterally
Others	Cross seat belt, Side support, Head support, Scoliosis support, Groin abductor, Foldable footrest, Removable rain cover



34 cm

50 kg

A solid red silhouette of a world map is centered in the background of the page.

“  
**Manual Wheeled  
Bath Chair**



Chassis	Steel
Front Wheel	5" filler
Rear Wheel	5" filler
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Removable, PU
Footrest	Removable, Side Opening and Folding
Fabric	black leather
Backrest	May come out
Others	Calf Band and Toilet Seat
Model	Suitable for 18" 20" and 22" chassis model construction.



50-55 cm




40-46 cm




120-135 kg








16-18 kg



Different Rim Options



Chassis	Aluminum
Front Wheel	10 cm filling
Rear Wheel	10 cm filling, lockable
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Back Foldable and Detachable, PU
Footrest	Removable, Side Opening, Folding and High Adjustable
Fabric	Black trim (flame retardant)
Backrest	May come out
Others	Calf Band
Model	Suitable for 18" 20" and 22" chassismodel construction.

				
50-55 cm	40-46 cm	120-135 kg	15-17 kg	Different Rim Options



Chassis	Steel, Folding
Front Wheel	8 "x1.25" padding
Rear Wheel	24" PU padding (Aluminum, composite)
Seating Width	18" (46 cm)
Carrying Capacity	100-120 kg
Armrest	Removable, PVC
Footrest	Removable, Folding and Highly Adjustable and Angle Adjustable
Fabric	Black Leather
Backrest	Angle adjustable
Others	Calf Support, Hand Brake, Head Support, Toilet Seat
Model	Suitable for 18" 20" and 22" chassis model construction.



40-45 cm



44-46 cm



100-120 kg



16-19 kg



Chassis	Steel, Folding
Front Wheel	8 "x2" padding
Rear Wheel	24" PU padding (Aluminum, composite)
Seating Width	18" (46 cm)
Carrying Capacity	120-135 kg
Armrest	Fixed, PVC
Footrest	Çıkabilir, Katlanır ve Yüksek Ayarlı
Fabric	Black Leather
Backrest	Sabit
Others	Adjustable Angle, Calf Support, Hand Brake, Toilet Seat
Model	Suitable for 18" 20" and 22" chassismodel construction.



60-62 cm



45-50 cm



120-135 kg



15-17 kg





# Spare Parts



Pneumatic and Padded Front Wheels



Pneumatic and Padded Rear Wheels



Plastic Front-Rear Wheels



24 V DC Motors



24 V DC Joystick



Battery Types

## Spare Parts



Sponge Cushion



Foot Pallet



PU and PVC Type Armrest



Chargers



Engine Brushes



Brakes



# Medical Products



**Rolator**  
PW30



**Reverse Walker**  
PW40



**Walker**  
PW50



**Walker**  
PW60



**Wheeled Walker**  
PW70



**Walker**  
PM809



**Toilet Seat Riser**  
PB6500



**Toilet Seat Riser**  
PB6600



**Folding Walking Stick**  
PB70

MADE IN  
TURKEY



# Comfortable WHEELCHAIR



PM330



## Medical Products



**Armpit Crutch**  
PB10-S



**Armpit Crutch**  
PB11-M



**Armpit Crutch**  
PB12-L



**Kanedyen**  
PB60



**Walking Stick with Stool**  
PB7310



**OPiV**  
*Pirmax*



+49 174 2674267



+90 (533) 636 67 42



+90 (850) 308 67 42



[info@opiagroup.eu](mailto:info@opiagroup.eu)



**Germany Headoffice / Hauptstelle:** Ikarusallee 15 30179 Hannover DEUTSCHLAND / GERMANY



**Germany Warehouse / Lager :** Chemnitzer Straße 11-13 30952 Hannover-Ronnenberg/Empelde DEUTSCHLAND / GERMANY



**Turkey Headoffice:** 15 Temmuz Mah. 148004 Nolu Cad No:7 ŞEHİTKAMİL / GAZİANTEP



SCAN ME

